

freedom

Stimulators™

stim
wave



Programmer Wearable Antenna Assembly Freedom-8A/-4A Receiver



Freedom System offers the same therapy as other systems... without an implanted battery, the system is only 5% the size, and is the Most Minimally Invasive.



Don't lose out on therapy options. Program for Low Frequency, High Frequency*, Surge, DRG, or PNS.



Upgrade your externals to evolve with your lifestyle and cutting edge technology. No surgery needed for upgrades.

*CAUTION - Investigational device. Limited by Federal (or United States) law to investigational use. Pulse Rate Capabilities above 1500 Hz to 10,000 Hz are limited by Federal (or United States) to Investigational Use.

Indications for Use. The Freedom System is intended as the sole mitigating agent, or as an adjunct to other modes of therapy used in a multidisciplinary approach for chronic, intractable pain of the trunk and/or lower limbs; all peripheral nerves excluding the cranial nerves. The Freedom System is contraindicated for patients who are unable to operate the Freedom System, are poor surgical risks, are pregnant, have an occupational exposure to high non-ionizing radiation, are implanted with cardiac system, or diathermy. Refer to the Instructions for Use provided with the Freedom System or Stimwave.com for potential adverse effects, warnings, and precautions prior to using this product. Rx Only. Caution: Federal law (USA) restricts these devices to sale by or on the order of a physician. Wireless Pain Relief® is a registered trademark of Stimwave. Stimwave LLC is a distributor of Stimwave Technologies Inc. and StimQ LLC.


Stimwave LLC
1310 Park Central Blvd S
Pompano Beach, Florida 33064
(800) 965-5134
contact@stimwave.com
stimwave.com/pns



04-00512-2
MD0707-03-4
© Stimwave LLC 2017

TREAT YOUR ARM & WRIST PAIN

With the Freedom Stimulator System
Less Pain and More Living™



wireless pain relief®

FREEDOM STIMULATION SOLUTION




ULNAR & RADIAL PAIN

The ulnar nerve runs from the shoulder to the hand and is responsible for carrying nerve signals between the hand and the brain, enabling motion and feeling in the hand and forearm. Ulnar neuropathy is an inflammation or compression of the ulnar nerve, resulting in numbness, tingling, and pain. Estimates are that **40% of Americans experience some form of this neuropathy at some point in their lives.** Ulnar neuropathy can be caused by compressive injury, by leaning or resting on the elbows or a blow to the elbow or ulnar nerve entrapment. Symptoms include:

- Weakness/ Tenderness in Hand
- Sensitivity to Cold
- Tingling in Palm & 4th/5th Fingers
- Elbow Tenderness
- Burning Feeling in Hand, Arm or Fingers

You may have tried occupational therapy, non-steroidal anti-inflammatories, splints, steroid injections, or even non-entrapment surgery to restore function and alleviate pain. You may have tried a multitude of these treatment options that provided some short-term relief, only to find that the pain has returned. If none of these options have really provided long-term relief, peripheral nerve stimulation with the Freedom Stimulator System may be a viable option to let you start living again with less pain!



Why Stimulation?

Stimulation uses electrical signals to interrupt pain signals to the brain. A stimulator, the size of a very small wire with electrodes on it, is placed near the painful nerve to deliver electrical pulses. The intensity and frequency of the electrical current are able to alleviate your pain and help you get back to a normal life! You control the therapy as the Freedom System integrates Wireless Pain Relief® into your clothes and daily activity, so you can have less pain and more living!

FULL BODY PAIN RELIEF®

